



CIGS Lab System









多元スパッタ装置 Sputter coater

CIGS 蒸着装置 CIGS coater

CBD 成膜装置 CBD coater

I-V 電気出力測定装置 I-V Meas. tester



Substrate Cleaning

Mo Contact Films by Sputtering

CIGS Coating by co-evaporation

Buffer Layer Coating by CBD

i-ZnO and TCO by Sputtering

Metallic Contact by Sputtering

Cell Scribing

I-V Measurement

PV-Labo1 Process Flow

